DE

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Seite 2

Infrarot-Sender, koppelbar

EN

Operation Instructions

Page 16

Infrared transmitter, combinable

FR

Mode d'emploi

Page 38

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Inhalt

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Contents

Important Safety Instructions	17
Hearing Safety	17
Battery Safety and Disposal	17
Pacemaker Safety	16
Transmitter Safety Warnings	18
System Overview	19
Installation Procedures	20
Determine Coverage Area	20
Optional: Additional »PRO IR-202« for Extended Coverage	21
Control Panel	22
Connect Power	23
Connecting an Audio Source	23
Adjust the Audio Level	24
Combined Mode	25
Troubleshooting	26
Mounting Options	26
Technical Specifications »PRO IR-202«	27
Maintenance and Care	28
Warranty	28
Disposal	28

Important Safety Instructions

Receiver Safety Warnings

Hearing Safety

Caution!

Many Humantechnik receivers are designed to amplify sounds to a high volume level, which could potentially cause hearing damage if used improperly. To protect your hearing and the hearing of others:

1. Make sure the volume is turned down before putting on the earphone or headphone - only then adjust the volume to a comfortable level.
2. Set the volume level at the minimum setting that you need to hear.
3. If you experience feedback (a squealing or howling sound), reduce the volume setting and move the microphone away from the earphone or headphone.
4. Do not allow children or other unauthorized persons access to this product.

Battery safety and disposal

Caution!

Many receivers are supplied with non-rechargeable alkaline batteries. Do not attempt to recharge non-rechargeable batteries; they may explode, release dangerous chemicals, cause burns, or cause other serious harm to the user or product.

Pacemaker safety

Caution!

If you have a pacemaker or other medical device, make sure that you are using the receiver in accordance with safety guidelines established by your physician or the pacemaker/medical device manufacturer.

Transmitter Safety Warnings

Warning!

To reduce the risk of fire or electric shock, do not expose the system to rain or moisture. Do not use this apparatus near water. The system should not be exposed to dripping or splashing, and objects filled with liquids such as beverages should not be placed on the transmitter or receivers. Clean only with a dry cloth.

Servicing or attempting to service this device will void the warranty!

- Refer servicing to qualified personnel. Servicing is required when the system has been damaged in any way: if liquid has been spilled or objects have fallen into the unit, if the unit has been exposed to moisture, if the unit does not operate normally, or if the unit has been dropped.
- Do not block any ventilation openings. Install in accordance with the manufacturer's instructions.
- Do not install near any heat sources such as radiators, heat sinks, stoves, or other apparatus that produces heat.
- Use only attachments / accessories as specified by the manufacturer.
- Unplug the transmitter during lightning storms or when unused for long periods of time.
- Be advised that different operating voltages require the use of different types of line cord and attachment plugs. Check the voltage in your area and use the correct type.
- Use only the power supply provided by Humantechnik. Other power supplies may have similar specifications, but may not be equivalent in emissions ratings, in-rush current, etc.
- Use of an unapproved power supply may leave the device partially or completely inoperable, and will void the warranty.
- Protect the power cord from being walked on or pinched, particularly at plugs, receptacles, and near the power jack on the transmitter.
- The mains plug or a switchable appliance coupler can be used to disconnect the device, so the disconnecting plug / coupler should remain easily accessible.

System Overview

The »PRO IR-202« two-channel infrared transmitter combines infrared modulator and emitter technology into a single mountable enclosure - which reduces operating costs, eliminates the need for rack space and eases set-up.

The »PRO IR-202« is ideal for high-quality audio programs such as music, theater and audio description. The »PRO IR-202« will accept any line level, balanced or unbalanced audio inputs.

Infrared receivers detect the transmission and convert the light signal back into audio signals. The 2.3/2.8 MHz operating frequencies minimize interference caused by high-efficiency lighting equipment.

A further »PRO IR-202« emitter can be connected for additional coverage area. Power and audio signal are shared between the emitters via the according connection cables. All necessary cables are included.

No FCC license or radio approval is required with this equipment. It can be used anywhere in the world.

Attention:

This equipment has been tested and found to comply with the according EC-regulations. The

limits stated therein are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area may cause interference in which case the user will be required to correct the interference at his own expense.

Please note:

Sunlight, a plasma monitor or television can degrade the audio quality of the »PRO IR-202« transmission.

For best performance, the transmitter should be positioned as far away as possible from any exposure to the sun, plasma monitor or television.

Installation Procedures

Determine Coverage Area

When using the »PRO IR-202« transmitter in single-channel mode with the PR 22+ receiver, the system coverage area will exceed 186 m². »PRO IR-202« units automatically adjust to give the best coverage area possible in both 1-channel and 2-channel modes.

Figure 1 illustrates the typical coverage pattern for the »PRO IR-202« This can be affected by direct/indirect sunlight, reflections on walls and room construction.

Reflections of the infrared light from walls, ceilings, and floors may change these patterns.


Important: Remember to point the transmitter towards the listening audience.

Remember: Most objects block infrared light. The transmitter cannot be concealed behind walls, glass, curtains, etc. These patterns are the direct radiation pattern. The infrared radiation does not drop to zero outside the illustrated patterns; it decreases. It still may be useable at a greater distance, depending on the receiver sensitivity and the reflective characteristics of the room.

The »PRO IR-202« can be operated with or without the faceplate, if desired. The faceplate can be removed and painted to match the room decor. Lightly sand the face-plate with 400 grit sandpaper, apply primer and two coats of the desired final color.

Do not paint the infrared lens behind the faceplate of the »PRO IR-202«!

Figure 1:
Receiver coverage area with »PRO IR-202« transmitter in single channel mode.

Additional »PRO IR-202« for Extended Coverage

For larger facilities, an additional emitter can be used in order to increase the total coverage area. The »PRO IR-202« is capable of driving another »PRO IR-202« in a combined mode as a secondary emitter.

The coverage area can be adjusted to the according requirements, by setting the direction of the secondary transmitter into the same direction, in sequence or in a different angle to the primary transmitter.

Figure 2:

Example: Two parallel placed emitters with overlapping illumination patterns to cover larger listening areas.


Figure 3:

Example: Extended coverage area with a second emitter pointing in the same direction



Figure 4:


Above view »PRO IR-202«

Control panel

Figure 5:

Rear view »PRO IR-202«

Connect Power

WARNING:

Mains voltage must not fall below 100V AC, or system performance will be greatly reduced. The device may shut down until a sufficient power level is restored.

1. Plug power supply output cord into the "Power IN" (DC IN) of the »PRO IR-202«.
2. Plug the power supply unit into the mains power outlet.

Connecting an Audio Source

Automatic Channel Detection/ Shutdown

The »PRO IR-202« offers two carrier frequency channels: Ch 1 (2.3 MHz) and Ch 2 (2.8 MHz). It automatically detects the presence of audio on the microphone and RCA audio inputs, and transmits on either/both channels when an audio signal is present on that channel. If audio is present on both the microphone and RCA jack for a given channel, the audio signals will be mixed.

When no audio signal is present on a channel for three minutes, the »PRO IR-202« will shut down that channel.

Operating the »PRO IR-202« in single-channel mode provides a 40 percent increase in range over two-channel mode.

Line Level Source

The »PRO IR-202« will accept line level, balanced or unbalanced audio inputs.

Microphone Source

Plug an electret microphone (e. g. Humantechnik EH1205-0) into the 3.5mm "mic in" jack.

To assign the microphone to Ch 1 (2.3 MHz) use the mic channel selector switch to select "2.3". To assign the microphone to Ch 2 (2.8 MHz), select "2.8".


Connecting an XLR-OUT source to RCA-IN

If required, a symmetric line-out port of an audio source device can be connected with an RCA-IN of the »PRO IR-202« by means of an adapter cable (adapter cable not included, please see the wiring diagram Figure 6).

Figure 6:

Wiring scheme:

Symmetric Line-Out port (XLR) to »PRO IR-202« RCA Line-In adapter cable

1= SHIELD, 2= RED (SIGNAL IN PHASE), 3=BLACK (SIGNAL)

Adjust the Audio Level

Line Level Source

After an audio signal has been connected to the transmitter, adjust the audio source gain until the »PRO IR-202« audio-level indicator for that channel blinks periodically.

Microphone Source

Speak into the microphone from a normal distance and at a natural voice level. Use the microphone gain adjustment on the »PRO IR-202« to set the gain to the point where the audio level indicator on the selected channel blinks periodically. If the audio level indicator is always on or on most of the time, the level is too high. If the audio level indicator is never on, the level is too low.

Combined mode

As previously described, a further »PRO IR-202« can be added as a secondary transmitter, in order to achieve a larger coverage area within larger facilities.

Connecting the two »PRO IR-202«:

Power supply:

Use the included power connection lead to connect the „DC-Out“ socket of the first »PRO IR-202« (primary transmitter) with the „DC-In“ socket of the second »PRO IR-202« (secondary transmitter).


Audio signal:

Use the included audio connection lead (3,5 mm stereo jack, on both ends) to connect the „sync out“ port of the primary transmitter with the „sync in“ port of the secondary transmitter.

After the two »PRO IR-202« have been connected with each other, the Line-In level indicator lights on both transmitters should light up simultaneously.

Figure 7:

Connecting two »PRO IR-202« for combined mode

Troubleshooting

The »PRO IR-202« Ch1 or Ch2 “Audio-level Indicator” is not lit:

- Make sure the »PRO IR-202« is plugged in.
- Make sure the audio input is connected properly. See page 23 - 24.
- Make sure an active audio signal with an adequate level is being sent to the »PRO IR-202« transmitter.
- In combined mode: Check the „sync“ connection („sync out“ port of primary transmitter with „sync in“ port of secondary transmitter).

No Sound Reproduction by Receivers:

- Make sure the receivers and transmitters are using the same frequency.
- If some of the receivers work but others don't, check for bad batteries or earphones.
- If none of the receivers work, check to see if the power and audio input cables are connected to the transmitter. Verify that the ch1 or ch2 audio-level indicators blink every now and then.
- Check to see if the transmitter is connected properly to the audio source. The audio-level indicators should flash on channels that have audio.
- Make certain the transmitter is not covered by objects which would block infrared light.

Sound Through the Receivers is Weak and Noisy:

- Try adjusting the audio input level (see page 24). If the signal sounds okay, you may need to reposition the transmitter beam or add an additional »PRO IR-202« secondary transmitter to the system.

Buzzing or Humming Noise in Sound System:

- Check for ground loops or noise on the input signal.
- Check if optical interferences are present (see page 19).

If these instructions do not address your problem or the issue persists, please call Humantechnik / AUDIOropa or your authorized dealer or representative.

Mounting Options

Use the included wall- and ceiling fixture set to mount the transmitter as desired. The set consists of connecting rods, a bracket and a ball-head clamp. The double-ended connecting rod (featuring ¼“ bolts at both ends) serves to mount two »PRO IR-202« on top of each other in combination. The »PRO IR-202« can also be mounted on any camera tripod with a 1/4” threaded connector.

Technical Specifications »PRO IR-202«

Dimensions, Weight	L 25,2 cm x B 7,9 cm x H 3,8 cm w/o faceplate or H 6,4 cm w/faceplate; 0,4 kg
Color	Anthracite
Power Supply	Power supply with IEC line cord, 100-240V AC input, 50-60Hz, / 24 V DC output at 800 mA
DC Power Input	2.5mm ID barrel connector, center positive
Sleep/Power Save Mode	Shuts off carrier when no audio is present for 3 minutes
Modulation	FM Wideband, ± 50 kHz deviation max, 50 μ S pre-emphasis
Carrier Frequencies	2.3 MHz (Ch1) and 2.8 MHz (Ch2). Default at power on = carriers off. Carriers are automatically enabled upon presence of audio.
Emitter IR Power	0.7 W
Audio inputs	RCA jack for Ch1 and Ch2 accept line level, balanced or unbalanced audio Line in balanced/unbalanced switch: selects balanced or unbalanced line level audio for RCA jacks
Microphone input	3.5 mm, stereo jack with signal and bias connected to tip, electret microphone compatible (4 VDC bias supply with 2.7k ohm series resistor)
Microphone gain adjust	Rotary potentiometer
Microphone channel switch	Selects microphone input to Ch1 or Ch2
Audio indicators	Red LED blinks at nominal audio level. One per channel
Signal-to-Noise Ratio	70 dB, (line input)
Frequency Response	95Hz to 17.6kHz, -3dB re 1kHz (line inputs) 125Hz to 17.0kHz, -3dB re 1kHz (microphone input)
Total Harmonic Distortion	<1% (1 kHz, nominal deviation, line or microphone input)
Operating Requirements	0-50°C
Mounting Kit	Wall or Ceiling mount bracket, omnidirectional ball-head clamp, connecting rods
Approvals	CE, WEEE, RoHS
Compatible Receivers	RCI-102, IP-112, PR-20+, PR-22+

Maintenance and Care

The »PRO IR-202« does not require any maintenance. If the unit becomes dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents. Do not set up the unit where it will be exposed to full sunlight for long periods. In addition it must be protected against excessive heat, moisture and severe mechanical shock.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Warranty

The »PRO IR-202« is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer. The warranty covers the repair of the product and returning it to you free of charge during the warranty period. It is essential that you send in the product in its original packaging, so do not throw the packaging away. The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer as well as a copy of the sales slip.

Always specify the product serial number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the re-cycling of electric and electronic units.

You protect health and environment by the correct disposal of this product. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

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Sarabec Ltd